

Agilent 6460	Triple Quadrupole LC/MS
System Type	LC-MS/MS (Triple Quadrupole)
Mass Range	m/z 5 – 3000
Resolution	0.7 Da (Auto), 0.5 Da (Manual)
Mass Accuracy	Up to ± 0.1 Da
Dynamic Range	$\sim 10^6$
Scan Speed	Up to 12,500 Da/s
MRM Capability	Up to 4000 transitions (Dynamic MRM)
Ion Sources	ESI / APCI / APPI / Multimode
Polarity Switching	~ 30 ms
Collision Cell	Hexapole (CID)
Software	MassHunter

